

# EXTRACTING FORCEPS



English Pattern

A



English Pattern

B



English Pattern,  
Anatomical Handle

C



American  
Pattern

D



Deep Grip Forceps  
With Hole Pattern

E



Hole Pattern

F



G



## EXTRACTING FORCEPS (English Pattern)

Щипцы экстракционные (английская форма)  
 Daviers (Forme Anglaise)  
 Zahnzangen (Englische Form)  
 提取强制 (英语模式)

Щипці екстракційні (англійська форма)  
 Forceps de Extracción (Forme Inglesa)  
 استخراج الملقط (نمط اللغة الإنجليزية)



**Pi1-01**  
Fig. 1  
upper laterals  
and canines



**Pi1-02**  
Fig. 7  
upper  
premolars



**Pi1-03**  
Fig. 13  
lower  
premolars



**Pi1-04**  
Fig. 17  
upper molars,  
right



**Pi1-05**  
Fig. 18  
upper molars,  
left



**Pi1-06**  
Fig. 22  
lower roots



**Pi1-07**  
Fig. 51A  
upper roots



**Pi1-08**  
Fig. 67  
upper third  
molars



**Pi1-09**  
Fig. 74N  
lower molar



**Pi1-10**  
Fig. 79  
lower third  
molars



**Pi1-11**  
Set of 10 Pieces  
Fig. 1  
Fig. 7  
Fig. 13  
Fig. 17  
Fig. 18  
Fig. 22  
Fig. 51A  
Fig. 67  
Fig. 74N  
Fig. 79



160mm  
Witzel Forceps  
Lower roots, Universal  
拔除上下颌牙根使用



165mm  
Witzel Forceps  
Lower roots, Universal  
拔除上下颌牙根使用



**Pi1-12**  
Witzel Forceps  
With Gold  
Set of 3 Pcs



165mm  
Witzel Forceps  
Upper roots, Universal  
拔除上下颌牙根使用

# EXTRACTING FORCEPS (For Children)

- Щипцы экстракционные детские
- Щипці екстракційні дитячі
- Daviers - forme anglaise
- Fórceps de extracción - forma inglesa
- Pinzeper per estrazione per bambini
- استخراج الملقط (للأطفال)

## BABY TOOTH EXTRACTION FORCEPS

**Pi1-13**  
English Pattern  
Klein Fig. 1  
Pedodontic,  
Upper Incisors



**Pi1-14**  
English Pattern  
Klein Fig. 2  
Pedodontic Upper Premolars



**Pi1-15**  
English Pattern  
Klein Fig. 3  
Pedodontic Upper  
Molars Universal



**Pi1-16**  
English Pattern  
Klein Fig. 4  
Pedodontic Upper root



**Pi1-17**  
English Pattern Klein Fig. 5  
Pedodontic Lower Incisors.



**Pi1-18**  
English Pattern  
Klein Fig. 6  
Pedodontic Lower Molars



**Pi1-20**  
Set of 7 Pieces



**Pi1-19**  
English Pattern  
Klein Fig. 7  
Pedodontic Lower Molars,



**Pi1-21**  
Comellas Forceps  
Anatomical Handle



**Pi1-22**  
Koehler  
Luxation forceps  
right & left molars




**Pi1-23**  
Fig. 00  
Ottavi inferiori  
Lower third molars

# EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Щипцы екстракційні (англійська форма)
- Daviers (Forme Anglaise)
- Forceps de Extracción (Forme Inglesa)
- Zahnzangen (Englische Form)
- استخراج الملقط (نمط اللغة الإنجليزية)
- 提取强制 (英语模式)

## DEEP GRIPING EXTRACTING FORCEPS SET

Pi1-24



- No Fracture
- No Tissue Damage
- Extraction With Easy



**Pi1-25**  
Fig. 34  
Upper Anterior



**Pi1-26**  
Fig. 35  
Upper Universal



**Pi1-27**  
Fig. 321  
Lower Universal



**Pi1-28**  
Fig. 36A  
Lower Premolars  
Deep-Gripping



**Pi1-29**  
Fig. 67  
Upper Third Molars



**Pi1-30**  
Fig. 79 T  
Lower Third Molars

**A-TRAUMATIC EXTRACTION FORCEPS**  
Improved tapered beaks for complete root contact.  
Cross serrations for a never slip grasp.  
Easily enter the sub gingival.

Pi1-31

## Traumatic Extracting Forceps Dual Grip Teeth Serration Set of 6 Pcs



**Pi1-32**  
Fig. 336  
lower premolars  
deep-gripping



**Pi1-33**  
Fig. 1  
Upper Anterior



**Pi1-34**  
Fig. 35  
Upper  
Universal



**Pi1-35**  
Fig. 36  
Lower Premolars  
Deep-Gripping




**Pi1-36**  
Fig. 67  
Upper Third  
Molars



**Pi1-37**  
Fig. 79A  
Apex Forceps Lower  
Universal Deep Root Pick

# EXTRACTING FORCEPS (American Pattern)

- Щипцы экстракционные (американская форма)
- Daviers - forme americane
- Fórceps de extracción - forma americana
- 提取强制 (美国模式)
- Щипці екстракційні (американська форма)
- Pinzaper per estrazione - forma americana
- استخراج الملقط (النمط الأمريكي)



**Pi1-38**  
Fig. W1  
Winter upper centrals  
and canines



**Pi1-39**  
Fig. 10S  
upper roots



**Pi1-40**  
Fig. 65  
upper incisors  
and roots



**Pi1-41**  
Fig. 69 Tomas  
upper and  
lower roots



**Pi1-42**  
Fig. 88R Nevius  
upper molars,  
right



**Pi1-43**  
Fig. 88L Nevius  
upper molars,  
left



**Pi1-44**  
Fig. 101  
Hull  
Deciduous teeth



**Pi1-45**  
Fig. 151  
Cryer  
lower teeth,  
for children



**Pi1-46**  
Fig. 210S  
upper third  
molars




**Pi1-47**  
Fig. 217  
lower molars



**Pi1-48**  
Fig. 222  
lower third molars



**Pi1-49**  
Set of 11 Pcs



**Pi1-50**  
Fig. W1  
Winter upper centrals  
and canines



**Pi1-51**  
Fig. 15  
Harris  
lower molars



**Pi1-52**  
Fig. 16  
lower molars



**Pi1-53**  
Fig. 17  
lower molars



**Pi1-54**  
Fig. 18R  
Harris  
upper molars,  
right



**Pi1-55**  
Fig. 23  
lower molars



**Pi1-56**  
Fig. 32 Parnly  
upper premolars  
and molars



**Pi1-57**  
Fig. 53 R  
upper molars,  
right



**Pi1-58**  
Fig. 150  
Cryer  
upper incisors,  
premolars,  
roots



**Pi1-59**  
Fig. 222  
lower third molars



**Pi1-60**  
Fig. 210H  
upper third molars



**Pi1-61**  
Set of 11 Pcs

- Щипцы экстракционные (американская форма)
- Daviers - forme américaine
- Fórceps de extracción - forma americana
- 提取强制 (美国模式)

- Щипці екстракційні (американська форма)
- Pinzere per estrazione - forma americana
- استخراج الملقط (النمط الأمريكي)

## 普通拔牙钳



型号	名称
Pi1-62	成人拔牙钳 150# (上颌切牙、双尖牙及牙根)
Pi1-63	成人拔牙钳 151# (下颌切牙、双尖牙及牙根)
Pi1-64	成人拔牙钳 210# (上颌第三磨牙)
Pi1-65	成人拔牙钳 222# (下颌第三磨牙)
Pi1-66	成人残根钳 (上颌)
Pi1-67	成人残根钳 (下颌)
Pi1-68	成人残根钳 (上、下颌万能)
Pi1-69	成人拔牙钳 1# (上颌切牙及牙根)
Pi1-70	成人拔牙钳 8# (上颌双尖牙)
Pi1-71	成人拔牙钳 17# (下颌磨牙)
Pi1-72	成人拔牙钳 18L# (上颌左侧磨牙)
Pi1-73	成人拔牙钳 18R# (上颌右侧磨牙)
Pi1-74	成人拔牙钳 44# (下颌切牙及牙根)
Pi1-75	成人拔牙钳 46# (下颌切牙)
Pi1-76	成人拔牙钳 48# (下颌牙根)
Pi1-77	成人拔牙钳 65# (上颌切牙及牙根)
Pi1-78	成人拔牙钳 69# (上下颌牙根)

## 微创拔牙钳



型号 名称  
Pi1-79 | 微创儿牙钳-上颌万用残根钳



型号 名称  
Pi1-80 | 微创儿牙钳-下颌万用残根钳



型号 名称  
Pi1-81 | 微创成人钳-上颌万用钳



型号 名称  
Pi1-82 | 微创成人钳-下颌万用钳



型号 名称  
Pi1-83 | 微创成人钳-上颌前牙



型号 名称  
Pi1-84 | 微创成人钳-下颌前牙



型号 名称  
Pi1-85 | 微创成人钳-上颌磨牙



型号 名称  
Pi1-86 | 微创成人钳-下颌磨牙



型号 名称  
Pi1-87 | 微创残根钳-下颌残根钳 (加长)



型号 名称  
Pi1-88 | 微创残根钳-上颌残根钳 (加长)

# EXTRACTING FORCEPS (English Pattern)

- Щипцы экстракционные (английская форма)
- Daviers (Forme Anglaise)
- Zahnzangen (Englische Form)
- 提取强制 (英语模式)

- Щипці екстракційні (англійська форма)
- Forceps de Extracción (Forme Inglesa)
- استخراج الملقط (نمط اللغة الإنجليزية)



**Pi1-89**  
Set of 37 Pcs

- |                                 |                 |   |
|---------------------------------|-----------------|---|
| 01 - Extracting Forceps         | Fig.1           | 1 |
| 02 - Extracting Forceps         | Fig.7           | 1 |
| 03 - Extracting Forceps         | Fig.13          | 1 |
| 04 - Extracting Forceps         | Fig.17          | 1 |
| 05 - Extracting Forceps         | Fig.18          | 1 |
| 06 - Extracting Forceps         | Fig.18A         | 1 |
| 07 - Extracting Forceps         | Fig.19          | 1 |
| 08 - Extracting Forceps         | Fig.20          | 1 |
| 09 - Extracting Forceps         | Fig.22          | 1 |
| 10 - Extracting Forceps         | Fig.29          | 1 |
| 11 - Extracting Forceps         | Fig.33          | 1 |
| 12 - Extracting Forceps         | Fig.51          | 1 |
| 13 - Extracting Forceps         | Fig.51A         | 1 |
| 14 - Extracting Forceps         | Fig.67A         | 1 |
| 15 - Extracting Forceps         | Fig.79          | 1 |
| 16 - Extracting Forceps         | Fig.136         | 1 |
| 17 - Extracting Forceps         | Fig.137         | 1 |
| 18 - Extracting Forceps         | Fig.53 L        | 1 |
| 19 - Extracting Forceps         | Fig.53 R        | 1 |
| 20 - Extracting Forceps         | Fig.65          | 1 |
| 21 - Root Elevators             | Fig.1B          | 1 |
| 22 - Root Elevators             | Fig.3B          | 1 |
| 23 - Root Elevators             | Fig.5B          | 1 |
| 24 - Root Elevators             | Fig.4B SA Apexo | 1 |
| 25 - Root Elevators             | Fig.5B SA Apexo | 1 |
| 26 - Cotton Tweezers            |                 | 1 |
| 27 - Mouth Mirrors With Handles | Fig.5           | 1 |
| 28 - Explorers                  |                 | 1 |
| 29 - Cement Spatulas            |                 | 1 |
| 30 - Filling Instruments        |                 | 1 |
| 31 - Filling Instruments        |                 | 1 |
| 32 - Filling Instruments        |                 | 1 |
| 33 - Filling Instruments        |                 | 1 |
| 34 - Gum Scissors               |                 | 1 |
| 35 - Excavators                 |                 | 1 |
| 36 - Bone Curettes              |                 | 1 |
| 37 - Dental Aspirated Syringes  |                 | 1 |



**Pi1-90** Available in Gold Coated Set of 25 Pcs  
**Pi1-91** Available in Polish Set of 25 Pcs

- |                                 |                 |   |
|---------------------------------|-----------------|---|
| 01 - Extracting Forceps         | Fig.13          | 1 |
| 02 - Extracting Forceps         | Fig.17          | 1 |
| 03 - Extracting Forceps         | Fig.18          | 1 |
| 04 - Extracting Forceps         | Fig.19          | 1 |
| 05 - Extracting Forceps         | Fig.20          | 1 |
| 06 - Extracting Forceps         | Fig.22          | 1 |
| 07 - Extracting Forceps         | Fig.29          | 1 |
| 08 - Extracting Forceps         | Fig.33          | 1 |
| 09 - Extracting Forceps         | Fig.51          | 1 |
| 10 - Extracting Forceps         | Fig.51A         | 1 |
| 11 - Extracting Forceps         | Fig.67A         | 1 |
| 12 - Extracting Forceps         | Fig.79          | 1 |
| 13 - Root Elevators             | Fig.1B          | 1 |
| 14 - Root Elevators             | Fig.3B          | 1 |
| 15 - Root Elevators             | Fig.4B SA Apexo | 1 |
| 16 - Root Elevators             | Fig.5B SA Apexo | 1 |
| 17 - Cotton Tweezers            |                 | 1 |
| 18 - Mouth Mirrors With Handles |                 | 1 |
| 19 - Explorers                  |                 | 1 |
| 20 - Cement Spatulas            |                 | 1 |
| 21 - Filling Instruments        |                 | 1 |
| 22 - Filling Instruments        |                 | 1 |
| 23 - Filling Instruments        |                 | 1 |
| 24 - Excavators                 |                 | 1 |
| 25 - Dental Aspirated Syringes  |                 | 1 |