



RETRACTORS - MOUTH GAGS

- Роторасширители
- Ecarteurs - Ouvre-bouches
- Vestibulumabhalter / Mundspanner
- 经销商-嘴塞
- Роторозширювачі
- Separadores - Abre-bocas
- الكاشطات - الكمامات الفم



Pi5-01 10.5cm
Pi5-02 13.5cm
Weitlaner
Fig. 1



Pi5-03 10.5cm
Pi5-04 13.5cm
Weitlaner
Fig. 2



Pi5-05
Mouth Gag Molt
9cm for (Children)



Pi5-06
Mouth Gag Molt
11.5cm for (Adults)



Pi5-07
Mouth Gag Molt
14cm for (Adults)



Pi5-08
Cheek Retractor
Long



Pi5-09
Cheek Retractor
Long



Pi5-10
14cm



Pi5-11
15cm

Cheek Retractor



Pi5-12
Unique Retractor
Length 17.5cm
angled 45°



Pi5-13
Unique Retractor
Length 17.5cm
angled 90°



Pi5-14
Cheek Retractor
Combination
21cm, 2 Pcs



Pi5-15
Sternberg
14 cm
Double Ended



Pi5-16
Sternberg
16 cm



Pi5-17
Sternberg
14 cm
Single Ended

RETRACTORS - MOUTH GAGS

- | | |
|----------------------------------|--------------------------|
| Роторасширители | Роторозширювачі |
| Ecarteurs - Ouvre-bouches | Separadores - Abre-bocas |
| Vestibulumabhalter / Mundspanner | الكاشطات - الكمامات الفم |
| 咬嘴器-嘴塞 | |



PI5-18
Stainless Steel Single Ended Cheek Retractors



PI5-19
Stainless Steel Double Ended Cheek Retractors



PI5-20
Cheek Retractor
Stainless Steel Long



PI5-21
Dental Lip Retractor



PI5-22
Left



PI5-23
Straight
Retractor Blue Titanium Coating



PI5-24
Right



PI5-25
Retractor
Blue Titanium Coating



PI5-26
Stainless Steel Double Black & Gold Ended Cheek Retractors

CRURBAN, University of Minnesota modified cheek retractor, is used to hold mucoperiosteal flaps, cheeks, lips and tongue away from the surgical area. 20.6mm



PI5-27
Urban Cheek Retractor

The Minnesota Retractor is used to hold mucoperiosteal flaps, cheeks, lips and tongue away from the surgical area. Retractor is 5.5" and 20.5mm wide and is made from surgical grade stainless steel for longevity and corrosion resistance.



PI5-28
University of Minnesota Retractor

The Camwood-Minnesota Cheek Retractor is used to hold mucoperiosteal flaps, cheeks, lips and tongue away from the surgical area. 5.9" long. Same design as the CRM, Minnesota Retractor with longer handle end. Made with surgical grade stainless steel.



PI5-29
Cawood Minnesota Cheek Retractor



PI5-30
Minnesota



PI5-31
Minnesota Cheek Retractor 15.5cm



PI5-32
Brunings Cheek Retractor 19cm




PI5-33
Cheek Retractor 20cm
Easy to Bend by Hand

RETRACTORS

Роторасширители
 Ecarteurs
 Wundhaken

Роторозширювачі
 Separadores
 الكاشطات



Pi5-43 30x11 mm
Pi5-44 30x16 mm
 Langenbeck Retractor
 21cm, hook



Middeldorpf Retractor
Pi5-45 17 cm
Pi5-46 22 cm



Pi5-47
 Farabeuf
 (120 X 12mm)
 Set of 2 Pcs



Pi5-48
 Seldin Retractor



Pi5-49
 Feldman
 Retractor
 14cm

RUBBER DAM CLAMPS FORCEPS

- Инструменты для системы коффердам
- Pincos porte-crampons a digue
- Kofferdam Klammerzangen
- 橡胶坝夹钳
- Инструменти для системи коффердам
- Pinzas de engrapar para dique de caucho
- السد المطاط الملتقط المشابك



Pi5-50

On request the Content of the rubber dam kit can be changed with other instruments.

Content:

- 01 Anisworth
- 02 Brewer
- 03 Clamps Stand
- 04 Rubber Dam Frame L
- 05 Rubber Dam Frame M
- 06 Rubber Dam Box Lower Case
- 07 Rubber Dam Box Upper Case
- 08 Rubber Dam Clamp Fig. 212
- 09 Rubber Dam Clamp Fig. 211
- 10 Rubber Dam Clamp Fig. 8
- 11 Rubber Dam Clamp Fig. 8A
- 12 Rubber Dam Clamp Fig. 14
- 13 Rubber Dam Clamp Fig. 14A
- 14 Rubber Dam Clamp Fig. 0
- 15 Rubber Dam Clamp Fig. 00
- 16 Rubber Dam Clamp Fig. 2
- 17 Rubber Dam Clamp Fig. 2A
- 18 Rubber Dam Clamp Fig. 7



Pi5-51
Fig. 15



Pi5-52
Fig. 18



Pi5-53
Fig. 24



Pi5-54
Fig. 25



Pi5-55
Fig. 26



Pi5-56
Fig. 28



Pi5-57
Fig. W3



Pi5-58
Fig. 4



Pi5-59
Fig. 56A



Pi5-60
Fig. 7



Pi5-61
Fig. 7A



Pi5-62
Fig. 8



Pi5-63
Fig. 210



Pi5-64
Fig. 211



Pi5-65
Fig. 212



Pi5-66
Fig. 6



Pi5-67
Fig. 9



Rubber Dam Kit Pi5-68

Anisworth 17cm, Ivory 17cm, Brewer 17.5cm
Rubber Dam Frame (Children), Ruber Dam Frame (Adult)
Ruber Dam Clamps 11 Pcs
Fig. 00, Fig.0, Fig.1, Fig.2A, Fig.3, Fig.7,
Fig.8, Fig.8A, Fig.9, Fig.212, Fig.14




Pi5-69
Anisworth



Pi5-70
Brewer



Pi5-71
Ivory



Pi5-73
Large Rubber Dam Frame

Pi5-72
Small Rubber Dam Frame



Pi5-74
Fig. 0



Pi5-75
Fig. 00



Pi5-76
Fig. 1



Pi5-77
Fig. 2A



Pi5-78
Fig. 3



Pi5-79
Fig. 7



Pi5-80
Fig. 8



Pi5-81
Fig. 8A



Pi5-82
Fig. 212



Pi5-83
Fig. 9



Pi5-84
Fig. 14



Pi5-85
Rubber Dam Punch



Pi5-86
Brewer



Pi5-87
Plier For Clamp




Pi5-88
Enc Clamp Kit 11 Pcs
Ruber Dam Clamps 9 Pcs
Fig.210, Fig.211, Fig.2, Fig.2A, Fig.207
Fig.201, Fig.202, Fig.203, Fig.204

RUBBER DAM CLAMPS

■ Резиновая дамба зажимы комплекты
■ Trousses de caoutchouc
■ Gummidammklammern
■ 橡胶夹

■ Гумова дамба затискачи комплекты
■ Conjuntos de abrazaderas de retención de caucho
■ المطاط المد المشابك مجموعات



All rubber dam clamps are made from surgical stainless steel. Matt finish of these Satin Steel Rubber Dam Clamps offset bright clinical lighting. Hand-crafted to absorb and dissipate light instead of reflecting it back. This section includes the Anterior Rubber Clamps.



Pi5-89
00

Flat jaws and high bow for small lower canines and first premolars

Pi5-90
6

Slightly offset jaws for labial caries on anterior teeth

Pi5-91
9

Offset flat jaws and stiff bow for anterior teeth

Pi5-92
9S


Offset flat serrated jaws for anterior teeth.

Pi5-93
212SA

Slightly offset staggered jaws for odd shaped anterior.

Pi5-94
212

Flat jaws for labial caries on anterior teeth.



All rubber dam clamps are made from surgical stainless steel. Matt finish of these Satin Steel Rubber Dam Clamps offset bright clinical lighting. Hand-crafted to absorb and dissipate light instead of reflecting it back. This section includes the Premolar Rubber Clamps.



Pi5-95
3

Flat narrow jaws for small lower molars.

Pi5-96
7

Flat jaws for lower molars

Pi5-97
7A

Jaws for large lower molars Slightly reversed festooned

Pi5-98
12A

Serrated offset jaws for lower right molars

Pi5-99
13A

Serrated offset jaws for lower left molars

Pi5-100
26

Slightly reversed festooned jaws for small lower molars

Pi5-101
28

Flat jaws for lower molars



Pi5-102
W3

Flat narrow jaws for small lower molars

Pi5-103
W7

Flat jaws for lower molars.

Rubber Dam Clamps -Premolar



Pi5-104
0

Flat jaws and high bow for longer premolars.

Pi5-105
1

Slightly festooned jaws for upper premolars.

Pi5-106
1A

Deeply festooned jaws for canines and premolars.

Pi5-107
2

Flat jaws for premolars.

Pi5-108
2A

Flat jaws for larger premolars.

Pi5-109
27N

Slightly festooned jaws for premolars.

Pi5-110
29

Flat jaws for lower premolars.

Pi5-111
209


Small, flat jaws with stiff bow for lower premolars

Pi5-112
W2

Flat jaws for premolars.

Rubber Dam Clamps -Molar (General)

All rubber dam clamps are made from surgical stainless steel. Matt finish of these Satin Steel Rubber Dam Clamps offset bright clinical lighting. Hand-crafted to absorb and dissipate light instead of reflecting it back. This section includes the Molar (General) Rubber Clamps.



Pi5-113
8A

Deeply festooned four point jaws for molars

Pi5-114
W8A

Slightly festooned jaws for molars.

Pi5-115
10

Flat jaws designed for small upper left and lower right molars.

Pi5-116
14

Flat jaws designed for small upper right & lower left molars

Pi5-117
14

Deeply festooned jaws for partially erupted molars

Pi5-118
14A

Deeply festooned jaws for partially erupted or irregularly shaped molars

Pi5-119
W14A

Deeply festooned jaws for partially erupted or irregularly shaped molars

Pi5-120
18

Flat jaws and stiff bow for molars

Pi5-121
24

Flat jaws for upper right and lower left molars

Pi5-122
25

Flat jaws for upper left and lower right molars

Pi5-123
27

Flat jaws for large molars

Pi5-124
56

Slightly festooned serrated jaws for large molars


Pi5-125
56S

Slightly festooned serrated jaws for large partially erupted molars

RUBBER DAM CLAMPS

■ Резиновая дамба зажимы комплекты
 ■ Trousses de caoutchouc
 ■ Gummidammklammern
 ■ 橡胶夹

■ Гумова дамба затискачи комплекти
 ■ Conjuntos de abrazaderas de retención de caucho
 ■ المطاط المد المشابك مجموعات



Rubber Dam Clamps -Upper Molar

All rubber dam clamps are made from surgical stainless steel. Matt finish of these Satin Steel Rubber Dam Clamps offset bright clinical lighting. Hand-crafted to absorb and dissipate light instead of reflecting it back. This section includes the Upper Molar Rubber Clamps.



Pi5-126

4
Flat narrow jaws for small upper molars



Pi5-127

W4
Flat narrow jaws for small upper molars



Pi5-128

5
Slightly festooned jaws for large upper molars.



Pi5-129

8
Deeply festooned jaws for upper molars



Pi5-130

26N
Flat jaws for upper molars.



Pi5-131

30
Offset jaws for upper right molars



Pi5-132

31
Offset jaws for upper left molars



Pi5-133

201
Slightly festooned jaws for upper molars



Pi5-134

205
Slightly festooned jaws for large upper molars.

Rubber Dam Clamps Organizing Board



Pi5-135
Rubber Dam Clamp Organizing Board, For 8 Clamps



Pi5-136
Rubber Dam Clamp Organizing Board, For 12 Clamps



Pi5-137
Stainless Steel Amalgam Well large



Pi5-138
Ivory without shields



Pi5-139
Tofflemire-Junior 45 mm



Pi5-140
Tofflemire-Universal 60 mm



Pi5-141
Matrix Ivory 8 6 mm



AMALGAM INSTRUMENTS TIP SIZES FOR DOUBLE-END

Pi5-142 - 1.5mm-2.0mm
Pi5-143 - 2.0mm-2.8mm
Pi5-144 - 2.0mm-3.2mm

RUBBER DAM CLAMPS

■ Резиновая дамба зажимы комплекты
■ Trousses de caoutchouc
■ Gummidammklammern
■ 橡胶夹

■ Гумова дамба затискачі комплекти
■ Conjuntos de abrazaderas de retención de caucho
■ المطاط السد المثابك مجموعات



Pi5-145 Fig. 15 **Pi5-146** Fig. 18 **Pi5-147** Fig. 24 **Pi5-148** Fig. 25 **Pi5-149** Fig. 26



Pi5-154 Fig. 138 **Pi5-155** Fig. 139 **Pi5-156** Fig. 200 **Pi5-157** Fig. 201 **Pi5-158** Fig. 202



Pi5-150 Fig. 28 **Pi5-151** Fig. 30 **Pi5-152** Fig. 31 **Pi5-153** Fig. 51



Pi5-159 Fig. 203 **Pi5-160** Fig. 204 **Pi5-161** Fig. 205 **Pi5-162** Fig. 3



Pi5-163 Fig. W3 **Pi5-164** Fig. 4 **Pi5-165** Fig. 56A **Pi5-166** Fig. 7 **Pi5-167** Fig. 7A



Pi5-172 Fig. 11 **Pi5-173** Fig. 12A **Pi5-174** Fig. 13A **Pi5-175** Fig. 14 **Pi5-176** Fig. 14A



Pi5-168 Fig. 8 **Pi5-169** Fig. 8A **Pi5-170** Fig. W8A **Pi5-171** Fig. 10



Pi5-177 Fig. 14AN **Pi5-178** Fig. 22 **Pi5-179** Fig. 27 **Pi5-180** Fig. 207



Pi5-181 Fig. 29 **Pi5-182** Fig. 206 **Pi5-183** Fig. 208 **Pi5-184** Fig. 209 **Pi5-185** Fig. 1



Pi5-190 Fig. 0 **Pi5-191** Fig. 00 **Pi5-192** Fig. 210 **Pi5-193** Fig. 211



Pi5-186 Fig. 1A **Pi5-187** Fig. 2 **Pi5-188** Fig. 2A **Pi5-189** Fig. W2A



Pi5-194 Fig. 212 **Pi5-195** Fig. 6 **Pi5-196** Fig. 9



Pi5-197 Fig. N06



Pi5-198 Fig. U67



Pi5-199 Fig. P1



Pi5-200 Fig. P2



Pi5-201 Fig. G-1



Pi5-202 Fig. G-2

CROWN REMOVER

Инструменты для удаления коронок
 Arrache-couronnes
 Teleskopkronenzangen
 皇冠去除

Инструменты для видалення коронок
 Arranca-corona
 مزيلات الساج

Fitting is possible on both ends



Pi5-203
Crown Remover,
For different level of intensity

Fitting is possible on both ends



Pi5-204
Crown Remover,
For different level of intensity

The inner spring can be remove, allow to control the stroke manually



Fig. 1 Fig. 2 Fig. 3 Fig. 4 Fig. 5 Fig. 6



Pi5-205
Crown Remover
Set of 3 Pcs



Pi5-206
Crown Remover
Set of 5 Pcs



Pi5-207
Crown Remover
Set of 6 Pcs



Pi5-208
Crown Remover
Set of 4 Pcs



Pi5-209
Offset Crown
Splitter
Remover
Elevator



Pi5-210
Christenson crown
remover
Straight, facial
use anterior



Pi5-211
Christenson crown
remover
Angled, for
posterior use



Pi5-212
Crown Spreaders



Pi5-213
Crown Removal Pliers
14cm



Rubber jaws
removes
Crown without
Damage
Crown Remover
Replacement
Soft Pads
Pi5-214

Free 10 Crown Remover Replacement Soft Pads


SCISSORS

- Ножницы десневые
- Ciseaux a gencives
- Zahnfleischscheren
- 口香糖剪刀
- Ножіці для ясен
- Tijeras para encias
- مقص صمغ



SCISSORS

- Ножницы десневые
- Ciseaux a gengives
- Zahnfleischscheren
- 口香糖剪刀
- Ножницы для ясен
- Tijeras para encias
- مقص صمغ



Иглодержатели
Porte-aiguilles
Nadelhalter

Голкотримачи
Porta-agujas
اصحاب ابرة

MICRO NEEDLE HOLDER/SCISSORS, FORCEPS AND SCALPEL HANDLES



Forceps
Pi5-256 14cm
Pi5-257 16cm
Pi5-258 18cm



Castroviejo Scissors TC
Straight
Pi5-259 14cm
Pi5-260 18cm



Castroviejo Scissors TC
Curved
Pi5-261 14cm
Pi5-262 18cm



Castroviejo Needle TC
Straight
Pi5-263 14cm
Pi5-264 18cm



Castroviejo Needle TC
Curved
Pi5-265 14cm
Pi5-266 18cm



Locks At 90°, 60°, 30°,
Straight, 120°, 150°, 180°,
Angles

BONE FILES



Pi5-267
Scalpel Handle
movable



Pi5-268
Scalpel Handle
Auto Push Blade



Pi5-269
Scalpel Handle
straight, #3



Pi5-270
Micro Surgical
Scalpel Handle



Pi5-271
Fig. 3
12.5cm



Pi5-272
Fig. 3
12.5cm



Pi5-273
Fig. 4
13.5cm



Pi5-274
Surgical Bone
Files-Miller
Bf45



Pi5-275
Surgical Bone
Files-Miller
Bf1x




Pi5-276
Surgical Bone
Files-Miller
Bf3x



Pi5-277
Surgical Bone
Files-Miller
Bf21


NEEDLE HOLDERS

Иглодержатели Голкотримачи
 Porte-aiguilles Porta-agujas
 Nadelhalter اصحاب ابرة



Pi5-281 10cm
 Pi5-282 12.5cm
 Pi5-283 14cm
 curved

Halstead Mosquito
 Pi5-278 10cm
 Pi5-279 12.5cm
 Pi5-280 14cm
 straight



Pi5-287 10cm
 Pi5-288 12.5cm
 Pi5-289 14cm
 Curved

Halstead Mosquito 1x2
 Pi5-284 10cm
 Pi5-285 12.5cm
 Pi5-286 14cm
 straight



Kelly-Rankin
 Pi5-290 16 cm
 straight
 Pi5-291 16 cm
 curved



Backhaus
 Pi5-292 11cm
 Pi5-293 13cm
 Pi5-294 14cm
 Pi5-295 15cm

T.C GOLD



Pi5-296 13cm
 Pi5-297 15cm
 Pi5-298 18cm
 Pi5-299 20cm

Crile-Ryder T.C



Pi5-300
 Halsey 12.5cm



Baumgartner



Pi5-309
 Derf 12 cm

Pi5-301 13 cm
 Pi5-302 14 cm
 Pi5-303 16 cm
 Pi5-304 18 cm
 Pi5-305 20 cm
 Pi5-306 25 cm
 Pi5-307 26 cm
 Pi5-308 30 cm



Pi5-310
 Delicate
 15 cm



Mayo-Hegar
 Pi5-311 13cm
 Pi5-312 16cm
 Pi5-313 18cm
 Pi5-314 20cm



Crile-Wood
 Pi5-315 15cm
 Pi5-316 17cm
 Pi5-317 20cm
 Pi5-318 22cm
 Pi5-319 25cm
 straight



Crile-Wood
 Pi5-320 15cm
 Pi5-321 17cm
 Pi5-322 20cm
 Pi5-323 22cm
 Pi5-324 25cm
 Curved



Pi5-325
 Mathieu
 14cm



Pi5-326
 Aruga
 16cm

BONE ROUNGER

DENTAL SYRINGES

Костные кусачки
Hohlmeißelzangen
Pincès-gouges

Кусачки для кісток
Pinzas-gubias
رونجيورس العظام

Шприцы стоматологические карпульные
Seringues à usage dentaire
Dentalspritzen
牙科注射器

Шприці стоматологічні карпульні
Jeringas dentales
المحقن الأسنان



Pi5-327
Surgical Bone Ronger - Friedman
14.5cm



Pi5-328
Surgical Bone Ronger - Mini Friedman
15° 12cm



Pi5-329
Surgical Bone Ronger - Micro Friedman
30° 14.5cm



Pi5-330
Surgical Bone Ronger - Micro Friedman
45° 14.5cm



Pi5-331
Surgical Bone Ronger - Micro Friedman
90° 14.5cm



Pi5-332
Ruskin Mini, Ronger Curved, 15cm, 5mm diameter tip



Pi5-333
Boehler, Ronger double action hinged, curved, 15cm, 2mm diameter tip



Pi5-334
Beyer, Ronger double action hinged, 18cm, 3mm diameter tip



Type C
Pi5-335 1.8ml
Pi5-336 2.2ml



3-Ring Non Folding
Pi5-340 1.8ml
Pi5-341 2.2ml



Pi5-346 1.8ml
Pi5-347 2.2ml
Angled



Pi5-337
Aspirated Dental Cartridge Syringe 1.8 ml



Aspirating Syringe
Pi5-342 1.8mm
Pi5-343 2.2mm



Aluminium - Dark
Pi5-348 1.8 ml
Pi5-349 2.2 ml



Pi5-338 1.8 ml
Pi5-339 2.2 ml
Aspirated Dental Cartridge Syringe Folding Style



Cartridge Syringe Folding
Pi5-344 1.8mm
Pi5-345 2.2mm



Aluminium - Blue
Pi5-350 1.8 ml
Pi5-351 2.2 ml



Projection System
Pi5-352 1.8ml
Pi5-353 2.2ml

Pi5-354
1001
Intraligamental Anesthesia Syringe in Black Colour

Pi5-355 1.8ml
Pistol-Ject
injects 0.20ml with every click Lever action for high injection pressure


Autoclavable



Pi5-356
Stieglitz Root Fragments



Pi5-357
Stieglitz Root Fragments



Pi5-358
Stieglitz Splinter Forceps
14cm 90°



Stieglitz Root Splinter Forceps



Pi5-359
Splinter Forceps Straight
18cm



Pi5-360
Splinter Forceps
18cm 45°