



# ORTHODONTIC INSTRUMENTS

■ Инструменты ортодонтические     ■ Инструменти ортодонтични  
■ Instrumentos Ortodonticos        ■ Instruments Orthodontiques  
■ Orthodontische Instrumente       ■ أدوات تقويم الأسنان  
■ 正畸仪器



**Pi6-01**  
Height Gauge 0.18  
Aluminum



**Pi6-02**  
Height Gauge 0.22  
Aluminum



**Pi6-03**  
Moveable  
Height Gauge



**Pi6-04**  
Bone Bracket Gauge



**Pi6-05**  
Endo Finger Ruler  
Span Measure Scale



**Pi6-06**  
Bone Thickness Gauge  
40mm, 18cm



**Pi6-07**  
Osteometer  
Measuring Forceps For Jaw  
Ridge Width measurement, scale  
on both sides, rustfree, autoclavable



**Pi6-08**  
Fig. 1  
Iwanson  
spring caliper  
For metal and  
porcelain.



**Pi6-09**  
Fig. 2  
Iwanson  
spring caliper  
For wax.



**Pi6-10**  
Castroviejo Calipers  
20mm, 45°  
9cm



**Pi6-11**  
Castroviejo Calipers  
20mm  
9cm



**Pi6-12**  
Generic 1Pc Dental Orthodontic Arch  
Forming Turret 0 Degree Without Torque  
available color



**Pi6-13**  
Direct Bond Bracket  
Holder Tweezers



**Pi6-14**  
Buccal tube tweezers  
(Extra Fine Tip)



**Pi6-15**  
Direct Bonding  
Bracket Holders  
13.5cm



**Pi6-16**  
Direct Bonding  
Bracket Holders  
13.6cm



**Pi6-17**  
Direct Bonding  
Bracket Holders  
13.5cm



**Pi6-18**  
Deluxe Bracket  
Holder

## ORTHODONTIC INSTRUMENTS

■ Инструменты ортодонтические      ■ Инструменти ортодонтични  
■ Instrumentos Ortodonticos      ■ Instruments Orthodontiques  
■ Orthodontische Instrumente      ■ ادوات تقويم الأسنان  
■ 正畸仪器

Available in Hollow Handle Ø9.5mm Handle



**Pi6-31**  
Wamkey Set of 3 Pcs Including Drill



**Pi6-27**  
Band seater  
auto cleavable  
Triangle Shape



**Pi6-28**  
Band seater  
auto cleavable  
Square Shape



**Pi6-29**  
Band seater  
auto cleavable  
Round Shape



**Pi6-30**  
Mershon Band  
Pusher  
14cm



**Pi6-32**



**Pi6-33**



**Pi6-34**

## PLIERS FOR ORTHODONTICS & PROSTHETICS

- |   |  |
|---|--|
| ■ Щипцы для ортодонтии и протезирования     | ■ Щипці для ортодонтії та протезування |
| ■ Pincers pour l'orthodontie et la prothèse | ■ Alicates para ortodoncia y prótesis  |
| ■ Zangen für Kieferorthopädie und Prothetik | ■ كمامشة لتقويم الأسنان بروتستيتيس     |
| ■ 牙齿矫正钳                                     |  |



**Pi6-35**  
Aderer  
Large wire bending  
and loop forming pliers  
13cm



**Pi6-36**  
Angle Plier  
Fig. 1  
13.5cm



**Pi6-37**  
Adams Plier  
Adams pliers for bending wire  
arches and bends right angles  
for wires up to .022" .025 tip  
to beak. Length 13cm



**Pi6-38**  
How Plier  
Universally usable tool for the  
multiband technology for wires  
up to .028". Length 14.5cm



**Pi6-39**  
How Plier curved  
Universally usable tool for  
the multiband technology for  
wires up to .028". Length 14.5cm



**Pi6-40**  
Weingart Plier medium  
Weingart plier that are ideal due to  
their curved tip for setting sheets and  
is capable of bending the ends.  
Length 14.5cm



**Pi6-41**  
Tweed Loop  
forming Pliers  
14cm



**Pi6-42**  
Nance Loop  
forming Pliers  
14cm



**Pi6-43**  
Band Removing Plier  
The curved blade of the band  
removing plier can strip easily  
and gently remove directly to  
the tooth. Length 14cm



**Pi6-44**  
Bracket Removing Plier curved for  
easy removal of brackets in the  
posterior region. Length 14.5cm



**Pi6-45**  
Bracket removing plier for  
comfortable and quick  
removal of the brackets.  
Length 13.5cm



**Pi6-46**  
Mini Ligature cutters T.C  
fine tip 13cm



**Pi6-47**  
Large  
Fluch-cut distal and cutting  
Pliers T.C with safety hold  
fitted with un-breakable,  
corrosion resistant and  
durable stainless steel tips 13cm



**Pi6-48**  
Small  
Fluch-cut distal and cutting  
Pliers T.C with safety hold  
fitted with un-breakable,  
corrosion resistant and  
durable stainless steel tips 13cm

## PLIERS FOR ORTHODONTICS & PROSTHETICS

- |   |                                      |
|---|--------------------------------------|
| Щипцы для ортодонтии и протезирования     | Щипці для ортодонтії та протезування |
| Pinces pour l'orthodontie et la prothèse  | Alicates para ortodoncia y prótesis  |
| Zangen für Kieferorthopädie und Prothetik | کماشة لتقويم الأسنان بروتستیس        |
| 牙齿矫正钳                                     |                                      |



**Pi6-49**  
Jarabak Pliers



**Pi6-50**  
Light Wire Pliers With Groove



**Pi6-51**  
Lingual Arch Forming Pliers



**Pi6-52**  
Ligatuer Bending plier



**Pi6-53**  
Bird Beak Pliers



**Pi6-54**  
Bracket Removing Plier  
Angled To - Inserts



**Pi6-55**  
Dental Orthodontic Pliers  
Crimpable Hook Placement



**Pi6-56**  
Crimpable Hook Placement



**Pi6-57**  
Rectangular Arch Forming Pliers



**Pi6-58**  
NiTi Three Jaw Pliers



**Pi6-59**  
Bird Beak Pliers With Cutter



**Pi6-60** Hammerhead  
NiTi Distal Cinch Back Pliers



**Pi6-61**  
Direct angled band plier  
anterior posterior bracket  
removing plier



**Pi6-62**  
Orthodontic Set of 21 Pcs



**Pi6-63**  
Orthodontic Set of 18 Pcs



**Pi6-64**  
Orthodontic Set of 10 Pcs

# PLIERS FOR ORTHODONTICS & PROSTHETICS

- Щипцы для ортодонтии и протезирования
 ■ Щипці для ортодонтії та протезування
- Pincers pour l'orthodontie et la prothèse
 ■ Alicates para ortodoncia y prótesis
- Zangen für Kieferorthopädie und Prothetik
 ■ کماشة لتقويم الأسنان بروستيثيس
- 牙齿矫正钳

## THERMO PLIERS & MATHIEU NEEDLE HOLDERS



**Pi6-65**  
Hole Punch Plier  
Creates half-moon cutouts in clear aligners to allow space for bonded buttons used for elastics



**Pi6-66**  
Vertical Punch Plier  
This versatile plier emphasizes on over-correction of rotations by creating an indentation at the mesial and/or lingual of a specific tooth. Can also be used to create a rotational couple using vertical indents on facial and lingual surfaces on opposite sides of a tooth.



**Pi6-67**  
Tear Drop Punch Plier  
Designed to create a notch at the gingival margin of clear aligners anywhere elastic hooks are required within the arch, making sure the elastic stays in place once the aligner is seated.



**Pi6-68**  
Horizontal Punch Plier  
Designed to produce an indentation to accent individual root torque. Can also be used to enhance retention of clear aligner or retainers by placing an indent into a tooth's undercut or below a bonded attachment.



**Pi6-69**  
Dimple Punch Plier  
Creates a tiny rounded dimple for spot contacts to accent individual tooth movements or added aligner retention. Typically used to accent movement on any tooth that requires a slight "nudge"



**Pi6-70**  
Eraser Plier  
The eraser plier is designed to remove and flatten any unwanted indentations placed on clear aligners by other punch pliers.



**Pi6-71**  
Groove thermal forming plier  
Designed for making grooves on clear appliances, to hold the elastic, prevent it falling off



**Pi6-72** ● Yellow  
Thermoplier  
Micro-Ramp Pliers  
Tip, torque or rotate smaller teeth laterally; Allow contact close to incisal edge  
On the correct location on the target tooth



**Pi6-73** ● Pink  
Thermoplier  
Small Circle Pliers  
Increase the size (depth) of small existing force point  
Location of the existing force point



**Pi6-74**  
Crown Beebe Scissors  
4.5 Str.



**Pi6-75**  
Mathieu Standard Tip



**Pi6-76**  
Elastomeric ligature placing forceps with groove 14cm



**Pi6-77**  
Mathieu Narrow Tip



**Pi6-78**  
Mathieu Smaha Ultra fine tip



**Pi6-81**  
Str.



**Pi6-82**  
Cvd.



**Pi6-79**  
Mathieu Standard Tip Tungsten Carbide



**Pi6-80**  
Mathieu Narrow Tip Tungsten Carbide

## PLIERS FOR ORTHODONTICS & PROSTHETICS

- |   |                                      |
|---|--------------------------------------|
| Щипцы для ортодонтии и протезирования     | Щипці для ортодонтії та протезування |
| Pinces pour l'orthodontie et la prothèse  | Alicates para ortodoncia y prótesis  |
| Zangen für Kieferorthopädie und Prothetik | كماشة لتقويم الأسنان بروتستيتيس      |
| 牙齿矫正钳                                     |                                      |



**Pi6-83**  
Dental Instruments  
Set of 24 Pcs

- 01 Minnesota Retractor
- 02 Mathieu Needle Holder
- 03 Tweezer
- 04 Probe
- 05 Scaler
- 06 Scaler
- 07 Scaler
- 08 Scaler
- 09 Scaler
- 10 Spatula
- 11 Mirror With Handle
- 12 Spatula
- 13 Scalpel Handles
- 14 Mead
- 15 Bone Ronger
- 16 Needle Holder
- 17 Crile Forcep
- 18 Towel Forcep
- 19 Iris Scissor
- 20 Goldman Scissor
- 21 Forcep
- 22 Adson Tweezer
- 23 Castroviejo Calipers
- 24 Amalgam Well



**Pi6-84**  
Dental Instruments  
Set of 25 Pcs

- 01 Scaler
- 02 Scaler
- 03 Scaler
- 04 Spatulas
- 05 Scaler
- 06 Mirror With Handle
- 07 Bone Files
- 08 Spatulas
- 09 Chiseles
- 10 Frazier
- 11 Minnesota
- 12 Dental Syringe
- 13 Castroviejo Calipers
- 14 Amalgam Well
- 15 Towel Forcep
- 16 Haemostatic Forcep
- 17 Needle Holder
- 18 Iris Scissor
- 19 Goldman Scissor
- 20 Scalpel Handle
- 21 B P Handle No.3
- 22 Adson Tweezer
- 23 College Tweezer
- 24 Bone Ronger